

**SIGVARIS TRADITIONAL MITTEN**  
Natural Rubber Quality



A-C Mitten

TRADITIONAL MITTEN	
<b>Models</b>	
A-C Mitten	⊙
<b>Colour</b>	
BEIGE	⊙

⊙ Series production

Packaging:  
1 pair in poly bag



Circumference in cm	small	medium	large
<b>cA</b>	16-19	19-22	22-25

Length in cm
1 Length

**Material composition:**  
70% polyamide  
30% elastodien

- Wash daily using a laundry bag
- Do not use fabric softener
- Do not place in dryer
- Do not iron
- Do not dry clean

Additional comments: A manifest lymphedema requires a complex, medically prescribed decongestif treatment. Regular examinations by a physician are necessary to maintain the result. In case of acute pain of skin consult your doctor immediately.

GANZONI & CIE. AG  
Gröblistrasse 8  
CH-9014 St.Gallen  
Telefon +41 (0)71 279 33 66  
Telefax +41 (0)71 274 29 29  
www.sigvaris-export.com

SIGVARIS, LIFE FOR LEGS and the leg icon are registered trademarks of Ganzoni & Cie AG, CH-9014 St.Gallen/Switzerland, in many countries worldwide.  
© 2010 Copyright by Ganzoni & Cie AG, St.Gallen/Switzerland

HZ/1000/03-2010/VZ/13982\_Ex

# SIGVARIS

## Arm Sleeves and Mitten

**SIGVARIS TRADITIONAL**  
Natural Rubber Quality

**SIGVARIS BASIC**  
Fine Quality

**SIGVARIS TRADITIONAL MITTEN**  
Natural Rubber Quality



## SIGVARIS TRADITIONAL

Natural Rubber Quality



- C-H without Mitten
- with support straps


- A-H with Mitten and thumb-piece
- with support straps

### INDICATIONS

- for recent lymphedemas that will subside rapidly and are soft
- post-mastectomy edema occurring late following radiotherapy
- post-operative and post-traumatic edema after completion of basic treatment
- polyarthritits in certain cases
- treatment by elastic compression for burns and extensive scarring

### CONTRAINDICATIONS

- unreduced lymphedema
  - inflammation
- No liability accepted in case of nonobservance of contraindications.



Circumference in cm	small	medium	large
cF	24-30	30-36	36-42
cD	20-24	24-28	28-32
cC	13-16	16-19	19-22

cF = Circumference of the arm at the thickest place  
cD = Circumference of the arm at mid forearm

Length in cm	
short	40-46 cm
long	46-52 cm



**Material composition:**  
70% polyamide  
30% elastodien

- Wash daily using a laundry bag
- Do not use fabric softener
- Do not place in dryer
- Do not iron
- Do not dry clean

TRADITIONAL				
CCL		2	3	4
<b>Models</b>				
C-H Arm Sleeve	without Mitten	●	■	■
A-H Arm Sleeve	with Mitten	●	■	■
<b>Colour</b>				
BEIGE		●		

CCL = Compression class  
● Series production and made-to-measure articles  
■ Made-to-measure articles

Packaging ½ pair per box including fitting/care instructions


## SIGVARIS BASIC

Fine Quality



- C-G without Mitten
- no support straps
- with and without grip top

- A-G with integral Mitten and thumb-piece
- no support straps
- with and without grip top



Circumference in cm	small	medium	large
cF	24-30	30-36	36-42
cD	20-24	24-28	28-32
cC	13-16	16-19	19-22

cF = Circumference of the arm at the thickest place  
cD = Circumference of the arm at mid forearm

BASIC				
CCL		1	2	
<b>Models</b>				
C-G Arm Sleeve	without Mitten	with soft top	⊙	⊙
		with grip top		⊙
A-G Arm Sleeve	with Mitten	with soft top	⊙	⊙
		with grip top		⊙
<b>Grip tops</b>				
Honey-comb grip top				⊙
<b>Colour</b>				
BEIGE				⊙

CCL = Compression class  
⊙ Series production

Packaging ½ pair per box including fitting/care instructions

Length in cm	
short	40-46 cm
long	46-52 cm



**Material composition:**  
70% polyamide  
30% elasthan

- Wash daily using a laundry bag
- Do not use fabric softener
- Do not place in dryer
- Do not iron
- Do not dry clean